

STRENGTH & SAFETY  
**IN EVERY DETAIL™**



 **Horton®**



# DRIVING HIGHER EXPECTATIONS

## OF STRENGTH, SAFETY, STYLE AND SERVICE.

2

There are no shortcuts or half measures here. For those who expect and require only the best, Horton redefines and reimagines safety, design and style with the ultimate customized ambulance. Our one-of-a-kind ambulances go beyond, employing application-specific and stress-tolerant materials, advanced engineering, more robust construction and the testing to prove it out.

Horton and its dedicated dealer network back this up with the kind of comprehensive warranties and service commitment that inspire confidence and make your decision easy.

Choosing a Horton ambulance is choosing to invest in your crew, your department and your mission, and to do it better every day.




### HOPS™ Featuring MBrace™

Our exclusive Horton Occupant Protection System™ (HOPS) featuring MBrace is the first ambulance patient compartment safety system that protects in both frontal impacts and rollover events. MBrace is a patent-pending, advanced frontal impact safety system that integrates an airbag into the best-in-class retractable harness for the attendant seat and side-facing positions, which include the CPR seat and squad bench. We also equip the attendant and CPR seats with the RollTek® Airbag System. In addition, progressive head foam is included for all patient compartment seating positions.

# THE HORTON DIFFERENCE



## Safety and Freedom to Work



## The Horton Legendary Body™

We've engineered a heavy-duty ambulance body with robust reinforcements—like a minimum of 31 fully boxed gussets per unit in the corners, walls and ceilings—keeping components connected and occupants protected. We also use more stress-tolerant aluminum for improved crack resistance with heavier gauge and wider extrusions for additional strength.




## Tested Beyond the Limits

We rigorously test strength, safety and integrity, employing third-party testing to be sure it's right. In fact, we exceed industry standards and requirements so that you can rest assured when you're on a run.



## Horton QuietGuard™

Our exclusive sound-deadening material is layered not just between frames, but over and inside them as well. We also add extra insulation by backing and lining our cabinets with felt and spray foam—creating a quieter environment for runs, so your crew can focus.



## Vi-Tech™ Mounting System

Exclusive Horton-engineered Vi-Tech body mounts employ neoprene inserts tested to 35Gs and connect our ambulance body outside the chassis frame rails for less vibration, better handling and more strength for smoother, quieter runs.



## HO<sub>2</sub>ton™ System

Our oxygen system employs a brass manifold to connect individual lines in the system and a uniquely robust regulator to meet the rigorous Horton standard of reliability. We build in a manual override, in case of power failure, and test every system at least two times the required pressure duration before installation.





# Horton® DISTINCTION IN EVERY DETAIL

## INNOVATIVE DESIGN YOU DIDN'T KNOW YOU NEEDED—AND CAN'T LIVE WITHOUT.

Your crews aren't the type to shy away from a challenge, and neither is Horton. While you're doing, we're always thinking. And when you're talking, we're listening. Energized by a challenge, our team is tasked with thinking beyond. Our passion for continual improvement and distinctive styling drives creativity and ingenuity inspired by feedback from the front lines. We bring cutting-edge advancement and design to your ambulance that further the safety and comfort of your people and those they serve.



## YOU KNOW WHAT YOU WANT.

We know we can deliver. Unique customization is standard operating procedure at Horton—that's why we make it easy. We're ready. Give us your challenge!






# THE HORTON DIFFERENCE



## Advanced Customization

Every ambulance design begins with a clean slate, and it's developed with your firsthand experience and our deep expertise. You're the expert on your department's needs, and we're the experts at getting it done.



## Design Lab

We foster a culture of discovery, and our engineers and designers are in constant pursuit of the next breakthrough—it's in our DNA. When you build an ambulance with Horton, you're not just choosing the latest and greatest—you're investing in your future.



## Distinctive Style


Our experts are ready to walk you through a design experience that will align the looks with your individualized needs. From flooring to cabinets to upholstery, we help you realize your unique vision that reflects the quality and esprit de corps of your department.

### THE HORTON DÉCOR BOOK


FOR MORE OPTIONS, OPEN THE TIP-OUT

At Horton, we make unique customization easy. See all the custom options at [hortonambulance.com](http://hortonambulance.com).


#### FLOORING FLECK



Moonstone




Sapphire Blue




Black Onyx




#### FLOORING PLATE



Gunmetal




Sapphire Blue




Brick




#### CARBON FIBER



Stainless Steel



Silver Stainless



Blue Stainless



## Cool-Tech®

Featuring the largest condenser in the industry, the powerful Cool-Tech system is available in three models. Cool-Tech intuitively monitors pressure inside the system and turns on only the number of fans needed to efficiently cool, alternating fans each time so they are used as evenly as possible. With a recessed pocket on top of the module, crews can enjoy additional body space for lighting, reduced noise and better protection from debris.



## Directional Safety Lighting

Horton's unique Directional Safety Lighting features a row of signal lights at the top of the ambulance that can be set to flash or illuminate from the center out, left to right or right to left. Operated from Power-Tech by Horton with one-touch activation and auto turn-off when the vehicle is taken out of park, Directional Safety Lighting eliminates a step when arriving at or leaving a scene.



## Custom Backlit Panels

Every Horton ambulance is ready to showcase your department with Custom Backlit Panels on the rear riser and side entry doorstep well that can illuminate a department name, crest or other messaging.




## Power-Tech by Horton™

Power-Tech by Horton is central command on every Horton ambulance, connecting to critical systems such as lighting, HVAC, the HO<sub>2</sub>rton System™ and Horton 360™, as well as vehicle diagnostics for components like Cool-Tech. The highly customizable system delivers flexibility and control features, with a programmable LCD touchscreen, to help save space on the front console and in the inhalation panel.




### THE HORTON DÉCOR BOOK

#### COMPOSITE CABINETS



White




Battleship  
Gray




Sandstone




#### COUNTERS



Midnight  
Pearl




Allspice  
Quartz




Gray  
Granite




#### LAMINATE



Brushed  
Aluminum




Brushed  
Copper




Aluminum  
Steel




#### UPHOLSTERY



Gunmetal



Dove



Port





## Horton 360™

Horton 360 delivers 360° views for video-based safety on the road and at the scene, providing an aerial view of the ambulance's surroundings, ensuring both crew and bystander safety. Horton 360 is accessible right from the driver's seat through Power-Tech by Horton, where crews can also customize inputs for additional ambulance video feeds.




## Custom Graphics

An ambulance is a point of pride for emergency departments and crews. Horton provides nearly unlimited options for customizing the appearance of your ambulance, with graphics ranging from individual paint schemes to full wraps and everything in between.




## Custom Rub-Rail Lighting

Add a distinctive flair to your Horton ambulance. Programmable rub-rail lights can increase visibility and distinction. They can be manually operated through Power-Tech by Horton or programmed to automatically operate in the overall scheme of the emergency lights.



Blue Ridge



Black

See all the custom options  
at [hortonambulance.com](http://hortonambulance.com).

Medium Duty

Type I

Type III




# DRIVEN TO SERVE

## BECAUSE WE KNOW YOU CAN'T SLOW DOWN

When you choose a high-performance ambulance, you expect high-performance coverage—and results. By partnering with some of the best dealers in the nation, we have a shared investment to provide unmatched training, assistance and solutions at the speed you need to do your job.

With Certified Emergency Vehicle Technicians on staff, our dealers receive Horton Training to continually improve their level of learning and service. One call to your Horton dealer takes care of it all, no matter the issue—even if it's a supplier part. Thanks to this ongoing collaboration, a Horton ambulance is a high-yield investment long after you hit the road.

Whatever the issue, your experienced Horton dealer will handle it. Whatever happens, we're all here to help.



**Our network of independent dealers,** backed by the Horton team and REV Finance, is ready as your first line of defense when you need support.



Scan to find your local dealer.

## HORTON COMPREHENSIVE COVERAGE

It only makes sense: A one-of-a-kind ambulance calls for one-of-a-kind support. Horton is there with warranties that deliver confidence and peace of mind.

## COMPETITIVE VEHICLE AND EQUIPMENT FINANCING

REV Financial Services offers affordable and flexible financing at competitive rates. Customers can choose from a variety of financing options and flexible payment schedules.

**Discover more at** [revgroup.com/rev-financial-services](https://revgroup.com/rev-financial-services).

## TO LEARN MORE OR BUILD YOUR OWN, VISIT [HORTONAMBULANCE.COM](https://hortonambulance.com)

Toll free	800.282.5113	3800 McDowell Road	<a href="mailto:info@hortonambulance.com">info@hortonambulance.com</a>
Direct	614.539.8181	Grove City, OH 43123	<a href="https://hortonambulance.com">hortonambulance.com</a>
Fax	614.539.8165		